

Endotoxin

Less than 1.0 EU per µg by the LAL method / rFC method.

Formulation

Supplied as 0.2 µm filtered solution in undefined.

Contact us for customized product form or formulation.

Shipping and Storage

This product is supplied and shipped on dry ice.

Please avoid repeated freeze-thaw cycles.

This product is stable after storage at:

- The product **MUST** be stored at -70°C or lower upon receipt;
- -70°C for 3 months under sterile conditions.

Discounts, Gifts,
and more!



www.acrobiosystems.com