

Biotinylated Mouse CD7 Protein, His Tag

Catalog # CD7-M82H3



BIOSYSTEMS
Acro

Surprise Inside!

Synonym

CD7, GP40, TP41, LEU-9, Tp40

Source

Biotinylated Mouse CD7 Protein, His Tag (CD7-M82H3) is expressed from human 293 cells (HEK293). It contains AA Gln 24 - Pro 150 (Accession # [P50283-1](#)). It is the Biotinylated form of Mouse CD7, His Tag (CD7-M52H3).

Predicted N-terminus: Gln 24

Molecular Characterization

CD7(Gln 24 - Pro 150)

P50283-1

Poly-his

This protein carries a polyhistidine tag at the C-terminus.

The protein has a calculated MW of 16.1 kDa. The protein migrates as 23-30 kDa when calibrated against [Star Ribbon Pre-stained Protein Marker](#) under reducing (R) condition (SDS-PAGE) due to glycosylation.

Labeling

The primary amines in the side chains of lysine residues and the N-terminus of the protein are conjugated with biotins using standard chemical labeling method. A standard biotin reagent (13.5 angstroms) is used in this product.

Protein Ratio

Passed as determined by the HABA assay / binding ELISA.

Purity

>90% as determined by SDS-PAGE.

Formulation

Lyophilized from 0.22 µm filtered solution in PBS, pH7.4 with trehalose as protectant.

Contact us for customized product form or formulation.

Reconstitution

Please see Certificate of Analysis for specific instructions.

For best performance, we strongly recommend you to follow the reconstitution protocol provided in the CoA.

Storage

For long term storage, the product should be stored at lyophilized state at -20°C or lower.

Please avoid repeated freeze-thaw cycles.


This product is stable after storage at:

- -20°C to -70°C for 12 months in lyophilized state;
- -70°C for 3 months under sterile conditions after reconstitution.

ACRO Quality Management System


- [QMS\(ISO, GMP\)](#)
- [Quality Advantages](#)
- [Quality Control Process](#)

SDS-PAGE



Biotinylated Mouse CD7 Protein, His Tag on SDS-PAGE under reducing (R) condition. The gel was stained with Coomassie Blue. The purity of the protein is greater than 90% (With [Star Ribbon Pre-stained Protein Marker](#)).

Bioactivity-ELISA



Immobilized Human SECTM1, His Tag (Cat. No. SE1-H5227) at 5 µg/mL (100 µL/well) can bind Biotinylated Mouse CD7 Protein, His Tag (Cat. No. CD7-M82H3) with a linear range of 0.001-0.625 µg/mL (QC tested).

Background

T-cell antigen CD7 (CD7) is also known as GP40, LEU-9, TP41 and Tp40. CD7 is a protein that in humans is encoded by the CD7 gene, this gene encodes a transmembrane protein which is a member of the immunoglobulin superfamily. CD7 has been shown to interact with PIK3R1. This protein is found on thymocytes and mature T cells. It plays an essential role in T-cell interactions and also in T-cell/B-cell interaction during early lymphoid development.

Discounts, Gifts,
and more!



www.acrobiosystems.com

