



ELISA Sample Dilution Buffer -2

Product information

Catalog Number: CEB-C305

Size: 15 mL/bottle

ELISA Sample Dilution Buffer-2 is a solution used to dilute and prepare test samples, facilitating target substance detection with ACRO ELISA reagents. It does not participate directly in the detection process.

Note: For the specific catalog number of kits compatible with this diluent, please contact Technical Support.

Storage

Keep the reagents stored at 2-8 °C. For opened kit and reconstituted reagents, can be stored for up to 30 days at 2-8 °C.

NOTE:

This dilution buffer product is for research use only and is not for usage in diagnostic or therapeutic applications.

Please do not use the reagents after the expiration date indicated on the label.

When ELISA Sample Diluent-2 is selected, choose one diluent exclusively for the entire experiment: use either the Sample Dilution Buffer provided with the kit, or ELISA Sample Diluent-2. Do not use both types, and do not switch diluents in the same assay

Manufactured and distributed by
ACRODiagnostics Inc.
TEL: 010-67855298-8147

Distributed by:

ACROBiosystems Inc.

US & Canada TEL: +1 800-810-0816

Asia & Pacific TEL: +86 400-682-252