

Synonym

Serum albumin,ALB,Alb

Source

Human Serum Albumin Protein, Flag Tag(HAS-H52C3) is expressed from human 293 cells (HEK293). It contains AA Asp 25 - Leu 609 (Accession # [P02768-1](#)).

Predicted N-terminus: Asp 25

Molecular Characterization



This protein carries a flag tag at the C-terminus.

The protein has a calculated MW of 68 kDa. The protein migrates as 65-70 kDa when calibrated against [Star Ribbon Pre-stained Protein Marker](#) under reducing (R) condition (SDS-PAGE) due to glycosylation.

Endotoxin

Less than 1.0 EU per µg by the LAL method / rFC method.

Purity

>90% as determined by SDS-PAGE.

Formulation

Lyophilized from 0.22 µm filtered solution in PBS, pH7.4 with trehalose as protectant.

Contact us for customized product form or formulation.

Reconstitution

Please see Certificate of Analysis for specific instructions.

For best performance, we strongly recommend you to follow the reconstitution protocol provided in the CoA.

Storage


For long term storage, the product should be stored at lyophilized state at -20°C or lower.

Please avoid repeated freeze-thaw cycles.

This product is stable after storage at:


- 20°C to -70°C for 12 months in lyophilized state;
- 70°C for 3 months under sterile conditions after reconstitution.

SDS-PAGE



Human Serum Albumin Protein, Flag Tag on SDS-PAGE under reducing (R) condition. The gel was stained with Coomassie Blue. The purity of the protein is greater than 90% (With [Star Ribbon Pre-stained Protein Marker](#)).

SEC-MALS




The purity of Human Serum Albumin Protein, Flag Tag (Cat. No. HAS-H52C3) is more than 85% and the molecular weight of this protein is around 60-90 kDa verified by SEC-MALS.

[Report](#)

Bioactivity-SPR

Discounts, Gifts,
and more!





Human Serum Albumin Protein, Flag Tag (Cat. No. HAS-H52C3) immobilized on CM5 Chip can bind Biotinylated Human FCGRT&B2M Heterodimer Protein, His,Avitag (Cat. No. FCM-H82W7) with an affinity constant of 202 nM as determined in a SPR assay (Biacore 8K) (QC tested).

Background

serum albumin (SA) is also known as ALB, which is the main protein of plasma and has a good binding capacity for water,Ca2+,Na+,K+,fatty acids,hormones, bilirubin and drugs.The main function of SA is the regulation of the colloidal osmotic pressure of blood. As Major zinc transporter in plasma, SA typically binds about 80% of all plasma zinc. A variant structure of albumin could lead to increased binding of zinc resulting in an asymptomatic augmentation of zinc concentration in the blood. Defects in serum albumin can cause familial dysalbuminemic hyperthyroxinemia which is a form of euthyroid hyperthyroxinemia that is due to increased affinity of serum albumin for T4. It is the most common cause of inherited euthyroid hyperthyroxinemia in Caucasian population.

