



Synonym

HLA-DRA1*0101 & HLA-DRB1*1101

Source

Biotinylated Human HLA-DRA1*01:01&HLA-DRB1*11:01 Monomer Protein(HL1-H82W6) is expressed from human 293 cells (HEK293). It contains AA Ile 26 - Glu 216 (HLA-DRA1*01:01) & Gly 30 - Lys 227 (HLA-DRB1*11:01) (Accession # [CAI2388006.1](#) (HLA-DRA1*01:01) & [CAI2388099.1](#) (HLA-DRB1*11:01)). Predicted N-terminus: Ile 26 & Gly 30

Molecular Characterization

Biotinylated Human HLA-DRA1*01:01&HLA-DRB1*11:01 Monomer Protein, produced by co-expression of HLA-DRA1*01:01 and HLA-DRB1*11:01, has a calculated MW of 29.2 kDa (HLA-DRA1*01:01) and 29.2 kDa (HLA-DRB1*11:01). Subunit HLA-DRA1*01:01 is fused with a polyhistidine tag at the C-terminus, followed by an Avi tag (Avitag™) and subunit HLA-DRB1*11:01 contains no tag. The protein migrates as 40-43 kDa and 35-38 kDa when calibrated against [Star Ribbon Pre-stained Protein Marker](#) under reducing (R) condition (SDS-PAGE) due to glycosylation.

Labeling

Biotinylation of this product is performed using Avitag™ technology. Briefly, the single lysine residue in the Avitag is enzymatically labeled with biotin.

Protein Ratio

Passed as determined by the HABA assay / binding ELISA.

Purity

>90% as determined by SDS-PAGE.
>90% as determined by SEC-MALS.

Formulation

Lyophilized from 0.22 µm filtered solution in PBS, pH7.4 with trehalose as protectant.
Contact us for customized product form or formulation.


Reconstitution

Please see Certificate of Analysis for specific instructions.
For best performance, we strongly recommend you to follow the reconstitution protocol provided in the CoA.

Storage


For long term storage, the product should be stored at lyophilized state at -20°C or lower.
Please avoid repeated freeze-thaw cycles.
This product is stable after storage at:
• -20°C to -70°C for 12 months in lyophilized state;
• -70°C for 3 months under sterile conditions after reconstitution.

SDS-PAGE



Biotinylated Human HLA-DRA1*01:01&HLA-DRB1*11:01 Monomer Protein on SDS-PAGE under reducing (R) condition. The gel was stained with Coomassie Blue. The purity of the protein is greater than 90% (With [Star Ribbon Pre-stained Protein Marker](#)).

SEC-MALS



The purity of Biotinylated Human HLA-DRA1*01:01&HLA-DRB1*11:01 Monomer Protein (Cat. No. HL1-H82W6) is more than 90% and the molecular weight of this protein is around 55-75 kDa verified by SEC-MALS.
[Report](#)




Biotinylated Human HLA-DRA1*01:01&HLA-DRB1*11:01 Monomer Protein (Peptide free, MALS verified)

Catalog # HL1-H82W6



Bioactivity-ELISA

Biotinylated Human HLA-DRA1*01:01&HLA-DRB1*11:01 Monomer Protein ELISA
0.1 µg of Biotinylated Human HLA-DRA1*01:01&HLA-DRB1*11:01 Monomer Protein per well



Immobilized Biotinylated Human HLA-DRA1*01:01&HLA-DRB1*11:01 Monomer Protein (Cat. No. HL1-H82W6) at 1 µg/mL (100 µL/well) on streptavidin (Cat. No. STN-N5116) precoated (0.5 µg/well) plate can bind Anti-HLA-DR Antibody, Human IgG1 with a linear range of 0.06-2 ng/mL (QC tested).

