

PE-Labeled Human HLA-A*02:01&B2M&p53 (HMTEVVRHC) Tetramer Protein

Catalog # HLP-HP2He



Synonym

HLA-A*0201 & B2M & p53 (HMTEVVRHC)

Source

PE-Labeled Human HLA-A*02:01&B2M&p53 (HMTEVVRHC) Tetramer Protein (HLP-HP2He) is expressed from human 293 cells (HEK293). It contains AA Ile 21 - Met 119 (B2M) & Gly 25 - Ile 308 (HLA-A*02:01) & HMTEVVRHC peptide (Accession # [P61769-1](#) (B2M) & [AAA59606.1](#) (HLA-A*02:01) & HMTEVVRHC).

Molecular Characterization

PE-Labeled Human HLA-A*02:01&B2M&p53 (HMTEVVRHC) Tetramer Protein is assembled by biotinylated monomer (HLP-H82Ef) and PE-labeled streptavidin.
Biotinylated Human HLA-A*02:01&B2M&p53 (HMTEVVRHC) Complex Protein is produced by co-expression of HLA and B2M loaded with p53 peptide. Biotinylated Human HLA-A*02:01&B2M&p53 (HMTEVVRHC) Complex Protein carries a polyhistidine tag at the C-terminus, followed by an Avi tag (Avitag™).

Conjugate

PE
Excitation Wavelength: 488 nm / 561 nm
Emission Wavelength: 575 nm

Formulation

Lyophilized from 0.22 µm filtered solution in PBS, 1% BSA, pH7.4 with trehalose as protectant.
Contact us for customized product form or formulation.

Reconstitution

Please see Certificate of Analysis for specific instructions.
For best performance, we strongly recommend you to follow the reconstitution protocol provided in the CoA.

Storage


For long term storage, the product should be stored at lyophilized state at -20°C or lower.
Please protect from light and avoid repeated freeze-thaw cycles.
This product is stable after storage at:

- -20°C to -70°C for 12 months in lyophilized state;
- -70°C for 3 months under sterile conditions after reconstitution.


ACRO Quality Management System

- [QMS\(ISO, GMP\)](#)
- [Quality Advantages](#)
- [Quality Control Process](#)

Bioactivity-ELISA




Immobilized PE-Labeled Human HLA-A*02:01&B2M&p53 (HMTEVVRHC) Tetramer Protein (Cat. No. HLP-HP2He) at 1 µg/mL (100 µL/well) can bind Anti-HLA class I Antibody, Human IgG1 (W6/32) with a linear range of 0.1-3 ng/mL (QC tested).



Immobilized PE-Labeled Human HLA-A*02:01&B2M&p53 (HMTEVVRHC) Tetramer Protein (Cat. No. HLP-HP2He) at 1 µg/mL (100 µL/well) can bind Anti-HLA-A2/p53R175H Antibody with a linear range of 0.005-0.625 µg/mL (QC tested).

PE-Labeled Human HLA-A*02:01&B2M&p53 (HMTEVVRHC) Tetramer Protein ELISA
0.1 µg of PE-Labeled Human HLA-A*02:01&B2M&p53 (HMTEVVRHC) Tetramer Protein per well



Immobilized PE-Labeled Human HLA-A*02:01&B2M&p53 (HMTEVVRHC) Tetramer Protein (Cat. No. HLP-HP2He) at 1 µg/mL (100 µL/well) can bind Anti-B2M Antibody, Human IgG1 with a linear range of 0.05-1 ng/mL (Routinely tested).

Background

TP53 has been recognized as a tumor suppressor. 50% of cancers carry a TP53 mutation while many others affect other pathway components. High-copy numbers of WT p53 peptide-MHC class I complexes were detected on tumor cells as compared to low copies on normal cells. The Human HLA-A*0201 p53 (HMTEVVRHC) complex protein is a complex of HLA-A*0201 of the MHC Class I, B2M, and HMTEVVRHC peptide of the p53.

Discounts, Gifts,
and more!



Acro



www.acrobiosystems.com