

Cynomolgus IL-31 Protein, His Tag, low endotoxin

Catalog # IL1-C5245



BIOSYSTEMS
Acro

Synonyms

IL-31, IL31, Interleukin-31

Source

Cynomolgus IL-31, His Tag (IL1-C5245) is expressed from human 293 cells (HEK293). It contains AA Leu 27 - Thr 163 (Accession # [G7PJ78-1](#)).

Predicted N-terminus: His

Molecular Characterization

Poly-his

IL-31(Leu 27 - Thr 163)
G7PJ78-1

This protein carries a polyhistidine tag at the N-terminus.

The protein has a calculated MW of 17.3 kDa. The protein migrates as 20-25 kDa when calibrated against [Star Ribbon Pre-stained Protein Marker](#) under reducing (R) condition (SDS-PAGE) due to glycosylation.

Endotoxin

Less than 0.01 EU per µg by the LAL method / rFC method.

Purity

>90% as determined by SDS-PAGE.

Formulation

Lyophilized from 0.22 µm filtered solution in PBS, pH7.4 with trehalose as protectant.

Contact us for customized product form or formulation.

Reconstitution

Please see Certificate of Analysis for specific instructions.

For best performance, we strongly recommend you to follow the reconstitution protocol provided in the CoA.

Shipping and Storage

This product is shipped at ambient temperature.

For long term storage, the product should be stored at lyophilized state at -20°C or lower.

Please avoid repeated freeze-thaw cycles.


This product is stable after storage at:

- -20°C to -70°C for 12 months in lyophilized state;
- -70°C for 3 months under sterile conditions after reconstitution.

ACRO Quality Management System


- [QMS\(ISO, GMP\)](#)
- [Quality Advantages](#)
- [Quality Control Process](#)

SDS-PAGE



Cynomolgus IL-31, His Tag on SDS-PAGE under reducing (R) condition. The gel was stained with Coomassie Blue. The purity of the protein is greater than 90% (With [Star Ribbon Pre-stained Protein Marker](#)).

Bioactivity-ELISA



Immobilized Cynomolgus IL-31, His Tag (Cat. No. IL1-C5245) at 1 µg/mL (100 µL/well) can bind Cynomolgus IL-31 RA Protein, Fc Tag (Cat. No. ILA-H525b) with a linear range of 5-78 ng/mL (QC tested).

Background

IL-31 is an inflammatory cytokine that helps trigger cell-mediated immunity against pathogens. Activates STAT3 and possibly STAT1 and STAT5 through the IL31 heterodimeric receptor composed of IL31RA and OSMR. It has also been identified as a major player in a number of chronic inflammatory diseases, including atopic dermatitis. May function in skin immunity. IL-31 is produced by a variety of cells, namely type 2 helper (TH2) T-cells. IL-31 sends signals through a receptor complex made of IL-31RA and oncostatin M receptor β (OSMR β) expressed in immune and epithelial cells. These signals activate three pathways: ERK1/2 MAP kinase, PI3K/AKT, and JAK1/2 signaling pathways.

Discounts, Gifts,
and more!



www.acrobiosystems.com



7/8/2026