

Biotinylated Human IL-18 R1 / CD218a Protein, Fc,Avitag™ (SPR & MALS verified)

Catalog # IL1-H82F9



Synonyms

IL18R1, IL1RRP, CD218a, IL18Ra

Source

Biotinylated Human IL-18 R1, Fc,Avitag (IL1-H82F9) is expressed from human 293 cells (HEK293). It contains AA Ala 19 - Arg 329 (Accession # [AAH69575.1](#)).

Predicted N-terminus: Ala 19

Molecular Characterization

IL-18 R1(Ala 19 - Arg 329) AAH69575.1	Fc(Pro 100 - Lys 330) P01857	Avi
--	---------------------------------	-----

This protein carries a human IgG1 Fc tag at the C-terminus, followed by an Avi tag (Avitag™).

The protein has a calculated MW of 63.9 kDa. The protein migrates as 70-90 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.

Labeling

Biotinylation of this product is performed using Avitag™ technology. Briefly, the single lysine residue in the Avitag is enzymatically labeled with biotin.

Protein Ratio

Passed as determined by the HABA assay / binding ELISA.

Purity

>95% as determined by SDS-PAGE.

Formulation

Lyophilized from 0.22 µm filtered solution in PBS, pH7.4 with trehalose as protectant.

Contact us for customized product form or formulation.

Reconstitution

Please see Certificate of Analysis for specific instructions.

For best performance, we strongly recommend you to follow the reconstitution protocol provided in the CoA.

Shipping and Storage

This product is shipped at ambient temperature.

For long term storage, the product should be stored at lyophilized state at -20°C or lower.

Please avoid repeated freeze-thaw cycles.


This product is stable after storage at:

- -20°C to -70°C for 12 months in lyophilized state;
- -70°C for 3 months under sterile conditions after reconstitution.

ACRO Quality Management System


- [QMS\(ISO, GMP\)](#).
- [Quality Advantages](#)
- [Quality Control Process](#)

SDS-PAGE




Biotinylated Human IL-18 R1, Fc,Avitag on SDS-PAGE under reducing (R) condition. The gel was stained with Coomassie Blue. The purity of the protein is greater than 95%.

SEC-MALS



The purity of Biotinylated Human IL-18 R1, Fc,Avitag (Cat. No. IL1-H82F9) is more than 85% and the molecular weight of this protein is around 145-175 kDa verified by SEC-MALS.

Bioactivity-SPR



Human IL-18, Tag Free immobilized on CM5 Chip can bind Biotinylated Human IL-18 R1, Fc,Avitag (Cat. No. IL1-H82F9) with an affinity constant of 0.628 μ M as determined in a SPR assay (Biacore 8K) (QC tested).

Background

Interleukin-18 receptor 1 (IL18R1) is also known as CD218 antigen-like family member A (CD218a), IL1 receptor-related protein (IL-1Rrp or IL1R-rp), CDw218a, IL18Ra, IL1RRP, which belongs to the interleukin-1 receptor family. IL18R1 contains three Ig-like C2-type (immunoglobulin-like) domains and one TIR domain. IL18R1 is receptor for interleukin 18 (IL-18). IL18R1 binds to the agonist leads to the activation of NF-kappa-B, but does not bind IL1A/interleukin-1 alpha or IL1B/interleukin-1 beta.



www.acrobiosystems.com