

Biotinylated Cynomolgus IL-33 Protein, His,Avitag™ (HPLC verified)

Catalog # IL3-C82E3



Synonym


IL33, DV27, C9ORF26, IL1F11, NFHEV, DKFZp586H0523, DVS27, NFEHEV, RP11-575C20.2

Source

Biotinylated Cynomolgus IL-33 Protein, His,Avitag (IL3-C82E3) is expressed from human 293 cells (HEK293). It contains AA His 109 - Ile 270 (Accession # [A0A2K5W3I1-1](#)).

Predicted N-terminus: His 109

Molecular Characterization



This protein carries a polyhistidine tag at the C-terminus, followed by an Avi tag (Avitag™).

The protein has a calculated MW of 21.9 kDa. The protein migrates as 27 kDa when calibrated against [Star Ribbon Pre-stained Protein Marker](#) under reducing (R) condition (SDS-PAGE) due to glycosylation.

Labeling

Biotinylation of this product is performed using Avitag™ technology. Briefly, the single lysine residue in the Avitag is enzymatically labeled with biotin.

Protein Ratio

Passed as determined by the HABA assay / binding ELISA.

Purity

>90% as determined by SDS-PAGE.

>95% as determined by SEC-HPLC.

Formulation

Lyophilized from 0.22 µm filtered solution in PBS, pH7.4 with trehalose as protectant.

Contact us for customized product form or formulation.

Reconstitution

Please see Certificate of Analysis for specific instructions.

For best performance, we strongly recommend you to follow the reconstitution protocol provided in the CoA.

Storage

For long term storage, the product should be stored at lyophilized state at -20°C or lower.

Please avoid repeated freeze-thaw cycles.


This product is stable after storage at:

- -20°C to -70°C for 12 months in lyophilized state;
- -70°C for 3 months under sterile conditions after reconstitution.

ACRO Quality Management System


- [QMS\(ISO, GMP\)](#)
- [Quality Advantages](#)
- [Quality Control Process](#)

SDS-PAGE




Biotinylated Cynomolgus IL-33 Protein, His,Avitag on SDS-PAGE under reducing (R) condition. The gel was stained with Coomassie Blue. The purity of the protein is greater than 90% (With [Star Ribbon Pre-stained Protein Marker](#)).

SEC-HPLC




The purity of Biotinylated Cynomolgus IL-33 Protein, His,Avitag (Cat. No. IL3-C82E3) was greater than 95% as determined by SEC-HPLC.

Bioactivity-ELISA




Immobilized Human IL-1RL1, Fc Tag (Cat. No. IL1-H5250) at 5 µg/mL (100 µL/well) can bind Biotinylated Cynomolgus IL-33 Protein, His,Avitag (Cat. No. IL3-C82E3) with a linear range of 0.2-2 ng/mL (QC tested).

Bioactivity-BLI



Loaded Biotinylated Cynomolgus IL-33 Protein, His,Avitag (Cat. No. IL3-C82E3) on SA Biosensor, can bind Human IL-1RL1, Fc Tag (Cat. No. IL1-H5250) with an affinity constant of 2.97 nM as determined in BLI assay (ForteBio Octet RED96e) (Routinely tested).



Loaded Human IL-1RL1, Fc Tag (Cat. No. IL1-H5250) on Protein A Biosensor, can bind Biotinylated Cynomolgus IL-33 Protein, His,Avitag (Cat. No. IL3-C82E3) with an affinity constant of 9.68 nM as determined in BLI assay (ForteBio Octet Red96e) (Routinely tested).

Background

Interleukin 33 (IL33) is known as C9orf26, DKFZp586H0523, DVS27, NF-HEV, NFEHEV, RP11-575C20.2, and is a cytokine belonging to the IL-1 superfamily. IL-33 induces helper T cells, mast cells, eosinophils and basophils to produce type 2 cytokines. IL-33 mediates its biological effects by interacting with the receptors

ST2 (aka IL1RL1) and IL-1 Receptor Accessory Protein (IL1RAP), activating intracellular molecules in the NF-κB and MAP kinase signaling pathways that drive production of type 2 cytokines (e.g. IL-5 and IL-13) from polarized Th2 cells. In vivo, IL-33 induces the expression of IL-4, IL-5, and IL-13 and leads to severe pathological changes in mucosal organs.



www.acrobiosystems.com