

Synonym

IL5Ra,CD125

Source

Mouse IL-5 R alpha Protein, His Tag(ILA-M52H7) is expressed from human 293 cells (HEK293). It contains AA Asp 18 - His 339 (Accession # [P21183](#)).
Predicted N-terminus: Asp 18

Molecular Characterization

IL-5 R alpha(Asp 18 - His 339)
P21183

Poly-his

This protein carries a polyhistidine tag at the C-terminus.

The protein has a calculated MW of 38.7 kDa. The protein migrates as 50-58 kDa and 75 kDa when calibrated against [Star Ribbon Pre-stained Protein Marker](#) under reducing (R) condition (SDS-PAGE) due to glycosylation.

Purity

>90% as determined by SDS-PAGE.

>90% as determined by SEC-MALS.

Formulation

Lyophilized from 0.22 µm filtered solution in PBS, pH7.4 with trehalose as protectant.

Contact us for customized product form or formulation.

Reconstitution

Please see Certificate of Analysis for specific instructions.

For best performance, we strongly recommend you to follow the reconstitution protocol provided in the CoA.

Storage


For long term storage, the product should be stored at lyophilized state at -20°C or lower.

Please avoid repeated freeze-thaw cycles.

This product is stable after storage at:


- 20°C to -70°C for 12 months in lyophilized state;
- 70°C for 3 months under sterile conditions after reconstitution.

SDS-PAGE



Mouse IL-5 R alpha Protein, His Tag on SDS-PAGE under reducing (R) condition. The gel was stained with Coomassie Blue. The purity of the protein is greater than 90% (With [Star Ribbon Pre-stained Protein Marker](#)).

SEC-MALS




The purity of Mouse IL-5 R alpha Protein, His Tag (Cat. No. ILA-M52H7) is more than 90% and the molecular weight of this protein is around 45-65 kDa verified by SEC-MALS.

[Report](#)

Bioactivity-SPR

Discounts, Gifts,
and more!





Mouse IL-5 R alpha Protein, His Tag (Cat. No. ILA-M52H7) immobilized on CM5 Chip can bind Mouse IL-5, His Tag (Cat. No. IL5-M52H3) with an affinity constant of 243 nM as determined in a SPR assay (Biacore 8K) (QC tested).

Background

Interleukin-5 receptor subunit alpha (IL5RA) is also known as IL-5 receptor subunit alpha, IL-5R subunit alpha and CD125, which is a single-pass type I membrane protein that belongs to the type I cytokine receptor family and type 5 subfamily. IL5R, a heterodimer of an alpha and a beta subunit, is expressed on eosinophils and basophils. IL5RA / CD125 is the receptor for interleukin-5 (IL5) and the alpha chain can bind to IL5. The WSXWS motif of IL5RA appears to be necessary for proper protein folding and thereby efficient intracellular transport and cell-surface receptor binding. The beta subunit of IL5R is common to the IL3, IL5 and GM-CSF receptors.

