

Synonym

IL-2 R beta & IL-2 R gamma,IL-2RB & IL-2RG

Source

Cynomolgus IL-2 R beta&IL-2 R gamma Heterodimer Protein, His Tag&Twin-Strep Tag(ILG-C5284) is expressed from human 293 cells (HEK293). It contains AA Ala 27 - Thr 240 & Leu 23 - Pro 255 (Accession # [Q38J85-1](#) & [A0A2K5VGJ6](#)).

Predicted N-terminus: Ala 27 & Leu 23

Molecular Characterization

IL-2 R beta (Ala 27 - Thr 240) Q38J85-1	Poly-his
IL-2 R gamma (Leu 23 - Pro 255) A0A2K5VGJ6	Twin-Strep

This protein carries a polyhistidine tag at the C-terminus, followed by a twin strep tag.

The protein has a calculated MW of 30.1 kDa & 34.1 kDa. The protein migrates as 38-43 kDa and 55-65 kDa when calibrated against [Star Ribbon Pre-stained Protein Marker](#) under reducing (R) condition (SDS-PAGE) due to glycosylation.

Endotoxin

Less than 1.0 EU per µg by the LAL method / rFC method.

Purity

>90% as determined by SDS-PAGE.

>90% as determined by SEC-MALS.

Formulation

Lyophilized from 0.22 µm filtered solution in PBS, pH7.4 with trehalose as protectant.

Contact us for customized product form or formulation.

Reconstitution

Please see Certificate of Analysis for specific instructions.

For best performance, we strongly recommend you to follow the reconstitution protocol provided in the CoA.

Storage


For long term storage, the product should be stored at lyophilized state at -20°C or lower.

Please avoid repeated freeze-thaw cycles.

This product is stable after storage at:

- 20°C to -70°C for 12 months in lyophilized state;
- 70°C for 3 months under sterile conditions after reconstitution.


SDS-PAGE



Cynomolgus IL-2 R beta&IL-2 R gamma Heterodimer Protein, His Tag&Twin-Strep Tag on SDS-PAGE under reducing (R) condition. The gel was stained with Coomassie Blue. The purity of the protein is greater than 90% (With [Star Ribbon Pre-stained Protein Marker](#)).

Bioactivity-ELISA

SEC-MALS




The purity of Cynomolgus IL-2 R beta&IL-2 R gamma Heterodimer Protein, His Tag&Twin-Strep Tag (Cat. No. ILG-C5284) is more than 90% and the molecular weight of this protein is around 90-120 kDa verified by SEC-MALS. [Report](#)



Cynomolgus IL-2 R beta&IL-2 R gamma Heterodimer Protein, His Tag&Twin-Strep Tag (MALS verified)

Catalog # ILG-C5284



Immobilized Cynomolgus IL-2, His Tag (Cat. No. IL2-C5249) at 5 µg/mL (100 µL/well) can bind Cynomolgus IL-2 R beta&IL-2 R gamma Heterodimer Protein, His Tag&Twin-Strep Tag (Cat. No. ILG-C5284) with a linear range of 1-78 ng/mL (QC tested).

Background

Both Interleukin-2 receptor subunit beta and Interleukin-2 receptor subunit gamma are receptor for interleukin-2. Common subunit for the receptors for a variety of interleukins. Interacts with SHB upon interleukin stimulation. Probably in association with IL15RA, involved in the stimulation of neutrophil phagocytosis by IL15. This beta subunit is involved in receptor mediated endocytosis and transduces the mitogenic signals of IL2. IL2R exists in 3 different forms: a high affinity dimer, an intermediate affinity monomer (beta subunit), and a low affinity monomer (alpha subunit). The high and intermediate affinity forms also associate with a gamma subunit.

Discounts, Gifts,
and more!

