

Monoclonal Anti-Inotuzumab Antibody, Mouse IgG2b (3F1)

Catalog # INB-MY2511



Source

Monoclonal Anti-Inotuzumab Antibody, Mouse IgG2b (3F1) is a recombinant mouse monoclonal antibody expressed in HEK293 cells.

Antibody Type

Recombinant Monoclonal

Clone

3F1

Isotype

Mouse IgG2b, Kappa

Host Species

Mouse

Specificity

Specifically recognizes Human Inotuzumab; This antibody is a non-neutralizing antibody.

Purification

Protein A purified / Protein G purified.

Aggregation

Less than 5%, as determined by SEC-MALS.

Concentration

Please refer to the Certificate of Analysis (CoA).

Form

Lyophilized

Formulation

Lyophilized from a 0.22 µm-filtered solution in PBS (pH 7.4), with trehalose as protectant.

Please contact us for customized product forms or formulations.

Reconstitution

Please refer to the Certificate of Analysis (CoA) for specific instructions.

Shipping

Lyophilized product is shipped at ambient temperature.

Storage

For long term storage, the product should be stored in a lyophilized state at -20°C or lower.

Please avoid repeated freeze-thaw cycles.

This product is stable after storage at:

- -20°C to -70°C for 12 months in lyophilized state;
- -70°C for 3 months under sterile conditions after reconstitution.


Notices

Product Specific Notices: For research use only.

ACRO Quality Management System

- [QMS\(ISO, GMP\)](#).
- [Quality Advantages](#)
- [Quality Control Process](#)

Bioactivity-ELISA



Immobilized Monoclonal Anti-Inotuzumab Antibody, Mouse IgG2b (3F1) (Cat. No. INB-MY2511) at 5 µg/mL, add Inotuzumab and then add Biotinylated Monoclonal Anti-Inotuzumab Antibody, Mouse IgG1 (7E5) (Cat. No. INB-BLMY2509) at 2 µg/mL. Detection was performed using HRP-conjugated Streptavidin (Acro, Cat. No. STN-NH913) (QC tested).

Background

Inotuzumab is a CD22-binding humanized monoclonal antibody. Inotuzumab is composed of the Fv (variable; antigen-binding) regions derived from the murine m5/44 anti-CD22 monoclonal antibody specific for the extracellular domain of human CD22 with human IgG4 heavy and kappa light chain constant (framework) regions.



www.acrobiosystems.com