

Mouse KLRG1 Protein, His Tag (MALS & SPR verified)

Catalog # KL1-M5249



Synonym

KLRG1

Source

Mouse KLRG1, His Tag (KL1-M5249) is expressed from human 293 cells (HEK293). It contains AA Gln 57 - Tyr 188 (Accession # [O88713-1](#)).

Molecular Characterization

Poly-his

KLRG1(Gln 57 - Tyr 188)
O88713-1

This protein carries a polyhistidine tag at the N-terminus.

The protein has a calculated MW of 17.1 kDa. The protein migrates as 24-27 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.

Endotoxin

Less than 1.0 EU per µg by the LAL method / rFC method.

Purity

>95% as determined by SDS-PAGE.

Formulation

Lyophilized from 0.22 µm filtered solution in PBS, pH7.4 with trehalose as protectant.

Contact us for customized product form or formulation.

Reconstitution

Please see Certificate of Analysis for specific instructions.

For best performance, we strongly recommend you to follow the reconstitution protocol provided in the CoA.

Storage

For long term storage, the product should be stored at lyophilized state at -20°C or lower.

Please avoid repeated freeze-thaw cycles.


This product is stable after storage at:

- -20°C to -70°C for 12 months in lyophilized state;
- -70°C for 3 months under sterile conditions after reconstitution.

ACRO Quality Management System


- [QMS\(ISO, GMP\)](#).
- [Quality Advantages](#)
- [Quality Control Process](#)

SDS-PAGE




Mouse KLRG1, His Tag on SDS-PAGE under reducing (R) condition. The gel was stained with Coomassie Blue. The purity of the protein is greater than 95%.

SEC-MALS




The purity of Mouse KLRG1, His Tag (Cat. No. KL1-M5249) is more than 85% and the molecular weight of this protein is around 19-27 kDa verified by SEC-MALS.

Bioactivity-ELISA



Immobilized Mouse KLRG1, His Tag (Cat. No. KL1-M5249) at 5 µg/mL (100 µL/well) can bind Human E-Cadherin (155-709) Protein, Fc Tag (Cat. No. ECD-H5256) with a linear range of 0.02-1.25 µg/mL (Routinely tested).

Bioactivity-SPR



Human E-Cadherin, Fc Tag (Cat. No. ECD-H5250) immobilized on CM5 Chip can bind Mouse KLRG1, His Tag (Cat. No. KL1-M5249) with an affinity constant of 0.159 µM as determined in a SPR assay (Biacore 8K) (QC tested).

Background

The co-inhibitory receptor killer-cell lectin like receptor G1 (KLRG1) is specifically expressed on NK cells and activated CD8+ T-cells and has been postulated to be a marker of senescence. KLRG1+ T cells are a major reason of chronic tissue damage in some autoimmune diseases such as systemic lupus erythematosus and rheumatoid arthritis. In tumors, tumor cells which express E-cadherin or N-cadherin bind to KLRG1 and inhibit the antitumor activity of T and NK cells. Thus, KLRG1 acts as an immuncheckpoint inhibitory receptor.

Discounts, Gifts,
and more!



www.acrobiosystems.com