

Human Klotho beta / KLB Protein, His Tag

Catalog # KLB-H52H4



BIOSYSTEMS
Acro

Synonym

betaKlotho, beta-klotho, BKL, KLB, klotho beta like, klotho beta-like protein

Source

Human Klotho beta Protein, His Tag (KLB-H52H4) is expressed from human 293 cells (HEK293). It contains AA Met 30 - Leu 997 (Accession # [Q86Z14-1](#)).

Predicted N-terminus: Met 30

Molecular Characterization

KLB(Met 30 - Leu 997) Q86Z14-1	Poly-his
-----------------------------------	----------

This protein carries a polyhistidine tag at the C-terminus.

The protein has a calculated MW of 112.9 kDa. The protein migrates as 120-130 kDa when calibrated against [Star Ribbon Pre-stained Protein Marker](#) under reducing (R) condition (SDS-PAGE) due to glycosylation.

Endotoxin

Less than 1.0 EU per µg by the LAL method / rFC method.

Purity

>90% as determined by SDS-PAGE.

Formulation

Lyophilized from 0.22 µm filtered solution in PBS, pH7.4 with trehalose as protectant.

Contact us for customized product form or formulation.

Reconstitution

Please see Certificate of Analysis for specific instructions.

For best performance, we strongly recommend you to follow the reconstitution protocol provided in the CoA.

Storage

For long term storage, the product should be stored at lyophilized state at -20°C or lower.

Please avoid repeated freeze-thaw cycles.


This product is stable after storage at:

- -20°C to -70°C for 12 months in lyophilized state;
- -70°C for 3 months under sterile conditions after reconstitution.

ACRO Quality Management System


- [QMS\(ISO, GMP\)](#).
- [Quality Advantages](#)
- [Quality Control Process](#)

SDS-PAGE




Human Klotho beta Protein, His Tag on SDS-PAGE under reducing (R) condition. The gel was stained with Coomassie Blue. The purity of the protein is greater than 90% (With [Star Ribbon Pre-stained Protein Marker](#)).

Bioactivity-ELISA




Immobilized Human Klotho beta Protein, His Tag (Cat. No. KLB-H52H4) at 10 µg/mL (100 µL/well) can bind Human FGF19, Fc Tag (Cat. No. FG9-H5253) with a linear range of 0.005-0.156 µg/mL (QC tested).



Immobilized Human Klotho beta Protein, His Tag (Cat. No. KLB-H52H4) at 5 µg/mL (100 µL/well) can bind Biotinylated Human FGF-21 Protein, His,Avitag (Cat. No. FG1-H82Q3) with a linear range of 0.02-0.156 µg/mL (Routinely tested).

Bioactivity-SPR



Human Klotho beta Protein, His Tag (Cat. No. KLB-H52H4) immobilized on CM5 Chip can bind Human FGF-21, His Tag (Cat. No. FG1-H5243) with an affinity constant of 67.1 nM as determined in a SPR assay (Biacore 8K) (Routinely tested).

Background

KLB (Klotho Beta) is a Protein Coding gene. Among its related pathways are RET signaling and HIV Life Cycle. GO annotations related to this gene include hydrolase activity, hydrolyzing O-glycosyl compounds and fibroblast growth factor binding. An important paralog of this gene is KL. Klotho Beta is a regulator in

multiple metabolic processes, while its role in cancer remains unclear. We found the expression of β Klotho was down-regulated in human hepatocellular carcinoma tissues compared with that in paired adjacent non-tumorous liver tissues. Hepatoma cells also showed decreased expression of β Klotho compared with normal hepatocyte cells. Reintroduction of β Klotho into hepatoma cells inhibited their proliferation.



www.acrobiosystems.com