

Cynomolgus NKp46 / NCR1 / CD335 Protein, His Tag (MALS verified)

Catalog # NC1-C52H4



Synonyms

NCR1, LY94, CD335, NK-p46, hNKp46

Source

Cynomolgus NKp46, His Tag (NC1-C52H4) is expressed from human 293 cells (HEK293). It contains AA Pro 22 - Arg 257 (Accession # [Q95JB9-1](#)).

Predicted N-terminus: Pro 22

Molecular Characterization

NKp46(Pro 22 - Arg 257) Q95JB9-1	Poly-his
-------------------------------------	----------

This protein carries a polyhistidine tag at the C-terminus.

The protein has a calculated MW of 28.6 kDa. The protein migrates as 34-40 kDa when calibrated against [Star Ribbon Pre-stained Protein Marker](#) under reducing (R) condition (SDS-PAGE) due to glycosylation.

Endotoxin

Less than 1.0 EU per µg by the LAL method / rFC method.

Purity

>90% as determined by SDS-PAGE.

>90% as determined by SEC-MALS.

Formulation

Lyophilized from 0.22 µm filtered solution in PBS, pH7.4 with trehalose as protectant.

Contact us for customized product form or formulation.

Reconstitution

Please see Certificate of Analysis for specific instructions.

For best performance, we strongly recommend you to follow the reconstitution protocol provided in the CoA.

Shipping and Storage

This product is shipped at ambient temperature.

For long term storage, the product should be stored at lyophilized state at -20°C or lower.

Please avoid repeated freeze-thaw cycles.


This product is stable after storage at:

- -20°C to -70°C for 12 months in lyophilized state;
- -70°C for 3 months under sterile conditions after reconstitution.

ACRO Quality Management System


- [QMS\(ISO, GMP\)](#).
- [Quality Advantages](#)
- [Quality Control Process](#)

SDS-PAGE




Cynomolgus NKp46, His Tag on SDS-PAGE under reducing (R) condition. The gel was stained with Coomassie Blue. The purity of the protein is greater than 90% (With [Star Ribbon Pre-stained Protein Marker](#)).

SEC-MALS



The purity of Cynomolgus NKp46, His Tag (Cat. No. NC1-C52H4) is more than 90% and the molecular weight of this protein is around 28-42 kDa verified by SEC-MALS.

Bioactivity-ELISA



Immobilized Cynomolgus NKp46, His Tag (Cat. No. NC1-C52H4) at 1 µg/mL (100 µL/well) can bind Monoclonal Anti-NKp46 Antibody, Human IgG1 with a linear range of 0.1-1 ng/mL (QC tested).

Background

Natural cytotoxicity triggering receptor 1 (NCR1) is also known as Natural killer cell p46-related protein (NK-p46), Lymphocyte antigen 94 homolog (LY94), CD antigen CD335, which belongs to the natural cytotoxicity receptor (NCR) family. NCR1 contains two Ig-like (immunoglobulin-like) domains. NCR1 interacts with CD247 and FCER1G. NCR1 / CD335 may contribute to the increased efficiency of activated natural killer (NK) cells to mediate tumor cell lysis.



www.acrobiosystems.com