

SARS-CoV-2 Nucleocapsid protein, His Tag (BA.3\*(S33V)/Omicron)

Catalog # NUN-C52Hv



Synonym

Nucleocapsid protein,NP,Protein N

Source

SARS-CoV-2 Nucleocapsid protein, His Tag (BA.3\*(S33V)/Omicron) (NUN-C52Hv) is expressed from human 293 cells (HEK293). It contains AA Met 1 - Ala 419 (Accession # [QHO62115.1](#) (P13L, GER 30-32del, S33V, R203K, G204R, S413R)). The mutations (P13L, GER 30-32del, S33V, R203K, G204R, S413R) were identified in the SARS-CoV-2 Omicron variant (Pango lineage: BA.3; GISAID clade: GRA; Nextstrain clade: 21M).  
Predicted N-terminus: Met 1

Molecular Characterization

This protein carries a polyhistidine tag at the C-terminus.  
The protein has a calculated MW of 47.1 kDa. The protein migrates as 50-60 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.

Endotoxin

Less than 1.0 EU per µg by the LAL method / rFC method.

Purity

>90% as determined by SDS-PAGE.

Formulation

Lyophilized from 0.22 µm filtered solution in PBS, 0.2 M Arginine, pH7.4 with trehalose as protectant.  
Contact us for customized product form or formulation.

Reconstitution

Please see Certificate of Analysis for specific instructions.  
*For best performance, we strongly recommend you to follow the reconstitution protocol provided in the CoA.*


Storage

For long term storage, the product should be stored at lyophilized state at -20°C or lower.  
*Please avoid repeated freeze-thaw cycles.*

This product is stable after storage at:

- -20°C to -70°C for 12 months in lyophilized state;
- -70°C for 3 months under sterile conditions after reconstitution.

SDS-PAGE



SARS-CoV-2 Nucleocapsid protein, His Tag (BA.3\*(S33V)/Omicron) on SDS-PAGE under reducing (R) condition. The gel was stained with Coomassie Blue. The purity of the protein is greater than 90%.

Bioactivity-ELISA




SARS-CoV-2 Nucleocapsid protein, His Tag (BA.3\*(S33V)/Omicron)

Catalog # NUN-C52Hv



SARS-CoV-2 Nucleocapsid protein, His Tag (BA.3\*(S33V)/Omicron) ELISA  
0.1 µg of SARS-CoV-2 Nucleocapsid protein, His Tag (BA.3\*(S33V)/Omicron) per well



Immobilized SARS-CoV-2 Nucleocapsid protein, His Tag (BA.3\*(S33V)/Omicron) (Cat. No. NUN-C52Hv) at 1 µg/mL (100 µL/well) can bind Anti-SARS-CoV-2 Nucleocapsid Antibody, Chimeric mAb, Human IgG1 (Cat. No. NUN-CH15) with a linear range of 0.5-2 ng/mL (QC tested).

Background

Nucleocapsid (N) protein is the most abundant protein found in coronavirus. CoV N protein is a highly immunogenic phosphoprotein important for viral genome replication and modulation of cell signaling pathways. It was first identified by a research team while they were screening for ADP-ribosylated proteins during coronavirus (CoV) infection (Grunewald M. E., et al. 2017, Virology; 517: 62-68). The array of diverse functional activities accommodated in N protein makes it more than a structural protein but also an interesting target in the development of antiviral therapeutics. Because of the conservation of N protein sequence and its strong immunogenicity, N protein of coronavirus is chosen as a diagnostic tool.

Discounts, Gifts,  
and more!

