

Cynomolgus PD-1 / PDCD1 Protein, His Tag, low endotoxin (MALS verified)

Catalog # PD1-C52H5



BIOSYSTEMS
Acro

Surprise Inside!

Synonym

PDCD1, PD1, CD279, SLEB2

Source

Cynomolgus PD-1, His Tag (PD1-C52H5) is expressed from human 293 cells (HEK293). It contains AA Leu 25 - Gln 167 (Accession # [B0LAJ3](#)).

Predicted N-terminus: Leu 25

Molecular Characterization

PD-1(Leu 25 - Gln 167)
B0LAJ3

Poly-his

This protein carries a polyhistidine tag at the C-terminus.

The protein has a calculated MW of 17.8 kDa. The protein migrates as 30-40 kDa when calibrated against [Star Ribbon Pre-stained Protein Marker](#) under reducing (R) condition (SDS-PAGE) due to glycosylation.

Endotoxin

Less than 0.01 EU per µg by the LAL method / rFC method.

Purity

>90% as determined by SDS-PAGE.

>95% as determined by SEC-MALS.

Formulation

Lyophilized from 0.22 µm filtered solution in PBS, pH7.4 with trehalose as protectant.

Contact us for customized product form or formulation.

Reconstitution

Please see Certificate of Analysis for specific instructions.

For best performance, we strongly recommend you to follow the reconstitution protocol provided in the CoA.

Storage

For long term storage, the product should be stored at lyophilized state at -20°C or lower.

Please avoid repeated freeze-thaw cycles.


This product is stable after storage at:

- -20°C to -70°C for 12 months in lyophilized state;
- -70°C for 3 months under sterile conditions after reconstitution.

ACRO Quality Management System


- [QMS\(ISO, GMP\)](#)
- [Quality Advantages](#)
- [Quality Control Process](#)

SDS-PAGE




Cynomolgus PD-1, His Tag on SDS-PAGE under reducing (R) condition. The gel was stained with Coomassie Blue. The purity of the protein is greater than 90% (With [Star Ribbon Pre-stained Protein Marker](#)).

SEC-MALS




The purity of Cynomolgus PD-1, His Tag (Cat. No. PD1-C52H5) is more than 95% and the molecular weight of this protein is around 25-35 kDa verified by SEC-MALS.

Bioactivity-ELISA




Immobilized Cynomolgus PD-1, His Tag (Cat. No. PD1-C52H5) at 10 µg/mL (100 µL/well) can bind Cynomolgus / Rhesus macaque PD-L1, Fc Tag (Cat. No. PD1-C5253) with a linear range of 0.039-1.25 µg/mL (QC tested).



Immobilized Cynomolgus PD-1, His Tag (Cat. No. PD1-C52H5) at 2 µg/mL (100 µL/well) can bind Nivolumab with a linear range of 0.1-1 ng/mL (Routinely tested).

Bioactivity-SPR



Biotinylated Cynomolgus / Rhesus macaque PD-L1 Protein, Fc, Avitag (Cat. No. PD1-C82F3) captured on Protein A Chip can bind Cynomolgus PD-1 Protein, His Tag (Cat. No. PD1-C52H5) with an affinity constant of 3.09 µM as determined in a SPR assay (Biacore 8K) (Routinely tested).

Background

Programmed cell death protein 1 (PD-1) is also known as CD279 and PDCD1, is a type I membrane protein and is a member of the extended CD28/CTLA-4 family of T cell regulators. PDCD1 is expressed on the surface of activated T cells, B cells, macrophages, myeloid cells and a subset of thymocytes. PD-1 has two ligands, PD-L1 and PD-L2, which are members of the B7 family. PD-L1 is expressed on almost all murine tumor cell lines, including PA1 myeloma, P815 mastocytoma, and

B16 melanoma upon treatment with IFN- γ . PD-L2 expression is more restricted and is expressed mainly by DCs and a few tumor lines. PD1 inhibits the T-cell proliferation and production of related cytokines including IL-1, IL-4, IL-10 and IFN- γ by suppressing the activation and transduction of PI3K/AKT pathway. In addition, coligation of PD1 inhibits BCR-mediated signal by dephosphorylating key signal transducer. In vitro, treatment of anti-CD3 stimulated T cells with PD-L1-Ig results in reduced T cell proliferation and IFN- γ secretion. Monoclonal antibodies targeting PD-1 that boost the immune system are being developed for the treatment of cancer.

Discounts, Gifts,
and more!



www.acrobiosystems.com



5/21/2026