

Mouse PD-L1 / B7-H1 Protein, Mouse IgG2a Fc Tag, low endotoxin (MALS verified)

Catalog # PD1-M52A2



BIOSYSTEMS
Acro

Synonym

PD-L1, CD274, B7-H1, PDCD1L1, PDCD1LG1

Source

Mouse PD-L1 Protein, Mouse IgG2a Fc Tag (PD1-M52A2) is expressed from human 293 cells (HEK293). It contains AA Phe 19 - Thr 238 (Accession # [NP_068693](#)).

Predicted N-terminus: Phe 19

Molecular Characterization

PD-L1(Phe 19 - Thr 238) NP_068693	mFc(Glu 98 - Lys 330) P01863
--------------------------------------	---------------------------------

[Other Tags and Version Biotin & Other Labeled Version](#)

This protein carries a mouse IgG2a Fc tag at the C-terminus.

The protein has a calculated MW of 51.8 kDa. The protein migrates as 65-85 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.

Endotoxin

Less than 0.01 EU per μg by the LAL method / rFC method.

Purity

>95% as determined by SDS-PAGE.

>95% as determined by SEC-MALS.

Formulation

Lyophilized from 0.22 μm filtered solution in PBS, pH7.4 with trehalose as protectant.

Contact us for customized product form or formulation.

Reconstitution

Please see Certificate of Analysis for specific instructions.

For best performance, we strongly recommend you to follow the reconstitution protocol provided in the CoA.

Storage

For long term storage, the product should be stored at lyophilized state at -20°C or lower.

Please avoid repeated freeze-thaw cycles.


This product is stable after storage at:

- -20°C to -70°C for 12 months in lyophilized state;
- -70°C for 3 months under sterile conditions after reconstitution.

ACRO Quality Management System


- [QMS\(ISO, GMP\)](#)
- [Quality Advantages](#)
- [Quality Control Process](#)

SDS-PAGE





Mouse PD-L1 Protein, Mouse IgG2a Fc Tag on SDS-PAGE under reducing (R) condition. The gel was stained with Coomassie Blue. The purity of the protein

SEC-MALS



The purity of Mouse PD-L1 Protein, Mouse IgG2a Fc Tag (Cat. No. PD1-M52A2) is more than 95% and the molecular weight of this protein is around

Bioactivity-ELISA



Immobilized Biotinylated Mouse PD-1 Protein, His,Avitag (Cat. No. PD1-M82E3) at 1 µg/mL (100 µL/well) on streptavidin (Cat. No. STN-N5116) precoated (0.5 µg/well) plate can bind Mouse PD-L1 Protein, Mouse IgG2a Fc Tag (Cat. No. PD1-M52A2) with a linear range of 0.06-2 ng/mL (QC tested).

Immobilized Mouse PD-L1 Protein, Mouse IgG2a Fc Tag (Cat. No. PD1-M52A2) at 5 µg/mL (100 µL/well) can bind Biotinylated Mouse PD-1, Fc,Avitag (Cat. No. PD1-M82F4) with a linear range of 0.313-5 µg/mL (Routinely tested).

Background

Programmed cell death 1 ligand 1 (PDL1) is also known as B7-H, B7H1, MGC142294, MGC142296, PD-L1, PDCD1L1 and PDCD1LG1, which is a member of the growing B7 family of immune molecules and is involved in the regulation of cellular and humoral immune responses. PDL1 is a cell surface immunoglobulin superfamily with two Ig-like domains within the extracellular region and a short cytoplasmic domain. This protein is broadly expressed in the majority of peripheral tissues as well as hematopoietic cells. Interaction between PDL1 and its receptors belonging to the CD28 family of molecules provide both stimulatory and inhibitory signals in regulating T cell activation and tolerance. PDL1 may inhibit ongoing T-cell responses by inducing apoptosis and arresting cell-cycle progression.

Discounts, Gifts,
and more!



www.acrobiosystems.com