

Biotinylated Monoclonal Anti-Sacituzumab Antibody, Mouse IgG1 (1H10)

Catalog # SAB-BLMY2500a



Source

Biotinylated Monoclonal Anti-Sacituzumab Antibody, Mouse IgG1 (1H10) is a recombinant mouse monoclonal antibody expressed in HEK293 cells.

Antibody Type

Recombinant Monoclonal

Clone

1H10

Isotype

Mouse IgG1, Kappa

Host Species

Mouse

Conjugate

Biotin

Specificity

Specifically recognizes Human Sacituzumab; This antibody is a neutralizing antibody.

Purification

Protein A purified / Protein G purified.

Aggregation

Less than 5%, as determined by SEC-MALS.

Concentration

Please refer to the Certificate of Analysis (CoA).

Form

Lyophilized

Formulation

Lyophilized from a 0.22 µm-filtered solution in PBS (pH 7.4), with trehalose as protectant.

Please contact us for customized product forms or formulations.

Reconstitution

Please refer to the Certificate of Analysis (CoA) for specific instructions.

Shipping

Lyophilized product is shipped at ambient temperature.

Storage

For long term storage, the product should be stored in a lyophilized state at -20°C or lower.

Please avoid repeated freeze-thaw cycles.

This product is stable after storage at:

- -20°C to -70°C for 12 months in lyophilized state;
- -70°C for 3 months under sterile conditions after reconstitution.


Notices

Product Specific Notices: For research use only.

ACRO Quality Management System

- [QMS\(ISO, GMP\)](#)
- [Quality Advantages](#)
- [Quality Control Process](#)

Bioactivity-ELISA



Immobilized Monoclonal Anti-Sacituzumab Antibody, Mouse IgG1 (1H7) (Cat. No. SAB-MY2501a) at 5 µg/mL, add Sacituzumab and then add Biotinylated Monoclonal Anti-Sacituzumab Antibody, Mouse IgG1 (1H10) (Cat. No. SAB-BLMY2500a) at 2 µg/mL. Detection was performed using HRP-conjugated Streptavidin (Acro, Cat. No. STN-NH913) (QC tested).

Background

Sacituzumab is a genetically engineered humanized monoclonal antibody directed against the human trophoblast cell-surface antigen 2 (TROP2, also known as TACSTD2) found on the surface of normal epithelial cells and a broad range of malignant solid tumor cells. The antibody is an IgG1 kappa immunoglobulin containing murine complementarity-determining region (CDR) sequences within human framework regions of the variable domains and human constant region sequences. Sacituzumab is composed of two heavy chains of 449 amino acids and two light chains of 214 amino acids.



www.acrobiosystems.com