

Human SCF (29-172) Protein, His Tag, low endotoxin (MALS verified)

Catalog # SCF-H52H4



Synonyms

KITLG, FPH2, KL-1, Kitl, MGF, SCF, SF, SHEP7, KL

Source

Human SCF (29-172) Protein, His Tag (SCF-H52H4) is expressed from human 293 cells (HEK293). It contains AA Cys 29 - Glu 172 (Accession # [AAH69797](#)).

Predicted N-terminus: Cys 29

Molecular Characterization

SCF(Cys 29 - Glu 172) AAH69797	Poly-his
-----------------------------------	----------

This protein carries a polyhistidine tag at the C-terminus.

The protein has a calculated MW of 18.2 kDa. The protein migrates as 20-30 kDa when calibrated against [Star Ribbon Pre-stained Protein Marker](#) under reducing (R) condition (SDS-PAGE) due to glycosylation.

The protein is designed as a dimer.

Endotoxin

Less than 0.01 EU per µg by the LAL method / rFC method.

Purity

>95% as determined by SDS-PAGE.

>90% as determined by SEC-MALS.

Formulation

Lyophilized from 0.22 µm filtered solution in PBS, pH7.4 with trehalose as protectant.

Contact us for customized product form or formulation.

Reconstitution

Please see Certificate of Analysis for specific instructions.

For best performance, we strongly recommend you to follow the reconstitution protocol provided in the CoA.

Shipping and Storage

This product is shipped at ambient temperature.

For long term storage, the product should be stored at lyophilized state at -20°C or lower.

Please avoid repeated freeze-thaw cycles.


This product is stable after storage at:

- -20°C to -70°C for 12 months in lyophilized state;
- -70°C for 3 months under sterile conditions after reconstitution.

ACRO Quality Management System


- [QMS\(ISO, GMP\)](#).
- [Quality Advantages](#)
- [Quality Control Process](#)

SDS-PAGE




Human SCF (29-172) Protein, His Tag on SDS-PAGE under reducing (R) condition. The gel was stained with Coomassie Blue. The purity of the protein is greater than 95% (With [Star Ribbon Pre-stained Protein Marker](#)).

SEC-MALS



The purity of Human SCF (29-172) Protein, His Tag (Cat. No. SCF-H52H4) is more than 90% and the molecular weight of this protein is around 45-60 kDa verified by SEC-MALS.

Bioactivity-ELISA



Immobilized Human CD117, Fc Tag (Cat. No. CD7-H5255) at 5 µg/mL (100 µL/well) can bind Human SCF (29-172) Protein, His Tag (Cat. No. SCF-H52H4) with a linear range of 0.02-5 ng/mL (QC tested).

Background

Stem Cell Factor is also known as SCF, kit-ligand, KL, steel factor, KITLG, FPH2, KL-1, Kitl, MGF, SCF, SF, or SHEP7, and is a cytokine that binds to the c-Kit receptor (CD117). SCF can exist both as a transmembrane protein and a soluble protein. This cytokine plays an important role in hematopoiesis (formation of blood cells), spermatogenesis, and melanogenesis. The soluble and transmembrane forms of the protein are formed by alternative splicing of the same R transcript. Soluble and transmembrane SCF is produced by fibroblasts and endothelial cells. Soluble SCF has a molecular weight of 18,5 KDa and forms a dimer. SCF plays an important role in the hematopoiesis during embryonic development. Sites where hematopoiesis takes place, such as the fetal liver and bone marrow, all express SCF. During development, the presence of the SCF also plays an important role in the localization of melanocytes, cells that produce melanin and control pigmentation. SCF plays a role in the regulation of HSCs in the stem cell niche in the bone marrow. SCF may be used along with other cytokines to culture HSCs and hematopoietic progenitors. The expansion of these cells ex-vivo (outside the body) would allow advances in bone-marrow transplantation, in which HSCs are transferred to a patient to re-establish blood formation.



www.acrobiosystems.com