

Mouse SOST / Sclerostin Protein, His Tag (HPLC verified)

Catalog # SOT-M5243



BIOSYSTEMS
Acro

Synonyms


SOST, VBCH

Source

Mouse SOST Protein, His Tag (SOT-M5243) is expressed from human 293 cells (HEK293). It contains AA Gln 24 - Tyr 211 (Accession # [Q99P68](#)).

Predicted N-terminus: His

Molecular Characterization



This protein carries a polyhistidine tag at the N-terminus.

The protein has a calculated MW of 22.5 kDa. The protein migrates as 30-35 kDa when calibrated against [Star Ribbon Pre-stained Protein Marker](#) under reducing (R) condition (SDS-PAGE) due to glycosylation.

Endotoxin

Less than 1.0 EU per µg by the LAL method / rFC method.

Purity

>90% as determined by SDS-PAGE.

>95% as determined by SEC-HPLC.

Formulation

Lyophilized from 0.22 µm filtered solution in PBS, 0.2 M Arginine, pH7.4 with trehalose as protectant.

Contact us for customized product form or formulation.

Reconstitution

Please see Certificate of Analysis for specific instructions.

For best performance, we strongly recommend you to follow the reconstitution protocol provided in the CoA.

Shipping and Storage

This product is shipped at ambient temperature.

For long term storage, the product should be stored at lyophilized state at -20°C or lower.

Please avoid repeated freeze-thaw cycles.


This product is stable after storage at:

- -20°C to -70°C for 12 months in lyophilized state;
- -70°C for 3 months under sterile conditions after reconstitution.

ACRO Quality Management System


- [QMS\(ISO, GMP\)](#).
- [Quality Advantages](#)
- [Quality Control Process](#)

SDS-PAGE




Mouse SOST Protein, His Tag on SDS-PAGE under reducing (R) condition. The gel was stained with Coomassie Blue. The purity of the protein is greater than 90% (With [Star Ribbon Pre-stained Protein Marker](#)).

SEC-HPLC




The purity of Mouse SOST Protein, His Tag (Cat. No. SOT-M5243) was greater than 95% as determined by SEC-HPLC.

Bioactivity-ELISA



Immobilized Mouse SOST Protein, His Tag (Cat. No. SOT-M5243) at 5 µg/mL (100 µL/well) can bind Human LRP-6 (20-630), Mouse IgG2a Fc Tag (Cat. No. LR6-H5253) with a linear range of 0.6-78 ng/mL (QC tested).



Immobilized Human LRP-6 (20-630), Mouse IgG2a Fc Tag (Cat. No. LR6-H5253) at 5 µg/mL (100 µL/well) can bind Mouse SOST Protein, His Tag (Cat. No. SOT-M5243) with a linear range of 0.005-0.156 µg/mL (Routinely tested).

Background

Sclerostin is a secreted glycoprotein with a C-terminal cysteine knot-like (CTCK) domain and sequence similarity to the DAN (differential screening-selected gene aberrative in neuroblastoma) family of bone morphogenetic protein (BMP) antagonists. Loss-of-function mutations in this gene are associated with an autosomal-recessive disorder, sclerosteosis, which causes progressive bone overgrowth. A deletion downstream of this gene, which causes reduced sclerostin expression, is associated with a milder form of the disorder called van Buchem disease. [provided by RefSeq, Jul 2008]



www.acrobiosystems.com